



Australia's premium manufacturer
of medical and laboratory products

Product Specification Sheet

Series	Algae growth chambers LED Lighting Temperature only
Model	agc-1620

1620 litre LED algae growth chamber
Temperature only



*Specifications subject to change without notice

Steridium Pty Ltd | Australia
www.steridium.com
sales@steridium.com
+61 7 55 931 352



Australia's premium manufacturer
of medical and laboratory products

Product Specification Sheet

Series Algae growth chambers
 LED Lighting
 Temperature only
Model agc-1620

1620 litre LED algae growth chamber
Temperature only

Chamber capacity	1620 Litres
Standard dual chamber design	1 set of shelves per door
Chamber width	705 mm
Chamber depth	600 mm
Chamber height	1810 mm
Single chamber design	Optional
Chamber width	1490 mm
Chamber depth	600 mm
Chamber height	1810 mm
Chamber material	Stainless Steel 304L
External width	1895 mm
External depth	762 mm
External height	2401 mm
External material	Durable epoxy powder-coated steel
Recommended clearance	75mm on sides and rear
Temperature range	4 to 50C
Temperature sensor	Pt100
Control accuracy	0.1 C
Temporal stability	0.2 C
Available shelf positions	38 (19 per door)
Shelves supplied	8 (4 per door)
Shelf dimensions	690 x 590 mm
Shelf spacing	90 mm
Shelf capacity	50 KG
Insulation type	Mineral Wool - asbestos free
Insulation thickness	50mm




Australia's premium manufacturer
of medical and laboratory products

Product Specification Sheet

Series Algae growth chambers
 LED Lighting
 Temperature only
Model agc-1620

1620 litre LED algae growth chamber
Temperature only

Compressor size	3/4 HP
Refrigerant	CFC-Free R134a
Lighting source	Bi-Directional (up/down) Waterproof LED Cassettes Lighting directly above and below entire shelf
Growth area	0.4m ² (per shelf)
Maximum intensity	200μmol/m ² /sec
Lighting colours	Lighting cassettes comprised of any combination of: Warm White Cool White Red Green Blue RGB Combination—Variable spectrum





Australia's premium manufacturer
of medical and laboratory products

Product Specification Sheet

Series Algae growth chambers
LED Lighting
Temperature only
Model agc-1620

1620 litre LED algae growth chamber
Temperature only

Electrical supply	Single phase, 220 - 240V, 50 - 60 Hz (alternatives available on request)
Element rating	2 kW
Circulating fan	2
Controller type	Multifunction intelligent digital PID-microprocessor
Readout	Digital LCD readout of set and actual temperature
Temperature alarms	High low deviation
Diurnal Control	Automatic switching between selected Day and Night Temperatures (and lights)
Country of origin	Australia



Australia's premium manufacturer
of medical and laboratory products

Product Specification Sheet

Series	Algae growth chambers LED Lighting Temperature only
Model	agc-1620

1620 litre LED algae growth chamber
Temperature only

Options / Accessories

- Single chamber design
- Additional shelves
- 8 Step programmable controller
- 500 Step programmable controller
- Inner viewing door - toughened glass
- Double glazed glass door
- Side access port - 50 mm
- Fan speed controller
- Internal power supply*
- RS232 or RS485 comms module
- Ethernet comms module
- ISPM 15 compliant crating for export
- Calibration certificate
- Co2 Control (ambient to 3000ppm)

*Maximum temperature limited to 50C