



Australia's premium manufacturer
of medical and laboratory products

Product Specification Sheet

Series Hot Air Sterilizers
Model h-2400

2400 litre hot air sterilizer





Australia's premium manufacturer
of medical and laboratory products

Product Specification Sheet

Series Hot Air Sterilizers
Model h-2400

2400 litre hot air sterilizer

| | |
|------------------------------|-----------------------------------|
| Chamber capacity | 2400 Litres |
| Standard dual chamber design | 1 set of shelves per door |
| Chamber width | 865 mm |
| Chamber depth | 900 mm |
| Chamber height | 1434 mm |
| Single chamber design | Optional |
| Chamber width | 1809 mm |
| Chamber depth | 900 mm |
| Chamber height | 1434 mm |
| Chamber material | Stainless Steel 304L |
| External width | 2214 mm |
| External depth | 1067 mm |
| External height | 2056 mm |
| External material | Durable epoxy powder-coated steel |
| Recommended clearance | 75mm on sides and rear |
| Temperature range | (Ambient + 7C) to 200C |
| Temperature sensor | Pt100 |
| Control accuracy | 0.1 C |
| Temporal stability | 0.2 C |
| Available shelf positions | 30 (15 per door) |
| Shelves supplied | 6 (3 per door) |
| Shelf dimensions | 855 x 900 mm |
| Shelf spacing | 87mm |
| Shelf capacity | 50 KG |
| Insulation type | Mineral Wool - asbestos free |
| Insulation thickness | 50mm |



Australia's premium manufacturer
of medical and laboratory products

Product Specification Sheet

Series Hot Air Sterilizers
Model h-2400

2400 litre hot air sterilizer

| | |
|--------------------|--|
| Circulating fan | 1 |
| Sterilizing cycle | 1 touch start to finish |
| Controller type | Multifunction intelligent digital PID-microprocessor |
| Readout | Digital LCD readout of set and actual temperature |
| Temperature alarms | High low deviation |
| Country of origin | Australia |

Options / Accessories

- Single chamber design
- Additional shelves
- RS232 or RS485 comms module
- Ethernet comms module
- ISPM 15 compliant crating for export
- Calibration certificate
- Swivel castors - Front lockable