



Australia's premium manufacturer  
of medical and laboratory products

## Product Specification Sheet

Series Hot Air Sterilizers  
Model h-30

30 litre hot air sterilizer





Australia's premium manufacturer  
of medical and laboratory products

## Product Specification Sheet

Series Hot Air Sterilizers  
Model h-30

30 litre hot air sterilizer

Chamber capacity	30 Litres
Chamber width	263 mm
Chamber depth	305 mm
Chamber height	360 mm
Chamber material	Stainless Steel 304L
External width	369 mm
External depth	505 mm
External height	618 mm
External material	Durable epoxy powder-coated steel
Recommended clearance	75mm on sides and rear
Temperature range	(Ambient + 7C) to 200C
Temperature sensor	Pt100
Control accuracy	0.1 C
Temporal stability	0.2 C
Available shelf positions	3
Shelves supplied	2
Shelf dimensions	258 x 280 mm
Shelf spacing	80mm
Shelf capacity	10 KG
Insulation type	Mineral Wool - asbestos free
Insulation thickness	50mm



Australia's premium manufacturer  
of medical and laboratory products

## Product Specification Sheet

Series Hot Air Sterilizers  
Model h-30

30 litre hot air sterilizer

Circulating fan	1
Sterilizing cycle	1 touch start to finish
Controller type	Multifunction intelligent digital PID-microprocessor
Readout	Digital LCD readout of set and actual temperature
Temperature alarms	High   low   deviation
Country of origin	Australia

### Options / Accessories

- Additional shelves
- RS232 or RS485 comms module
- Ethernet comms module
- ISPM 15 compliant crating for export
- Calibration certificate