



Australia's premium manufacturer  
of medical and laboratory products

## Product Specification Sheet

Series Peltier Cooled Waterbaths  
Model wp-40

40 litre peltier cooled waterbath



\*Optional stainless steel lid shown

\*Specifications subject to change without notice

Steridium Pty Ltd | Australia  
[www.steridium.com](http://www.steridium.com)  
[sales@steridium.com](mailto:sales@steridium.com)  
+61 7 55 931 352



Australia's premium manufacturer  
of medical and laboratory products

## Product Specification Sheet

Series Peltier Cooled Waterbaths  
Model wp-40

40 litre peltier cooled waterbath

Chamber capacity	40 Litres
Chamber length	623 mm
Chamber width	503 mm
Chamber depth	160 mm
Chamber material	Stainless Steel 304L
External length	898 mm
External width	554 mm
External depth	400 mm
External material	Durable epoxy powder-coated steel
Recommended clearance	50mm on sides and rear
Temperature range	(5C below Ambient) to 50C
Temperature sensor	Pt100
Control accuracy	0.1 C
Temporal stability	0.2 C
Insulation type	Mineral Wool - asbestos free
Insulation thickness	50mm
Electrical supply	Single phase, 220 - 240V, 50 Hz (alternatives available on request)
Element rating	1 kW
Countdown timer	99 hr 59 min - or continuous
Controller type	Multifunction intelligent digital PID-microprocessor
Readout	Digital LCD readout of set and actual temperature
Temperature alarms	High   low   deviation
Drain	Standard
Circulation pump	Standard
Country of origin	Australia



Australia's premium manufacturer  
of medical and laboratory products

## Product Specification Sheet

Series Peltier Cooled Waterbaths  
Model wp-40

40 litre peltier cooled waterbath

### Options / Accessories

- Stainless steel gabled lid
- 8 Step programmable controller
- ISPM 15 compliant crating for export
- Calibration certificate

\*Maximum temperature limited to 50C